

**Section:** Clinical

**Procedure No.** [#v1/MM/YYYY]

**PROCEDURE TITLE:** DNA extraction from mucosal tissue

**Review Officer:** Gastroenterology and  
Hepatology Consultant

**Review Summary:** v1

**Applicable To:** All Gastroenterology and  
Hepatology Department Laboratory Staff

**Date Introduced or Last Review Date:**  
11/2019

**Next Review Date:** 11/2022

**Authority:** Head of Clinical Function  
Testing, Gastroenterology and  
Hepatology

**Replaces:** New procedure

**Key Words:** DNA, mucosa

## PURPOSE

To ensure that the procedure is conducted safely and accurate results are obtained

## OUTCOME

The DNA extraction is carried out in a standardised and reproducible manner

## AUTHORISED TO UNDERTAKE THE PROCEDURE

Gastroenterology and Hepatology Laboratory Staff

## RISKS AND PRECAUTIONS

Appropriate PPE must be worn by laboratory technicians

## STEPS OF THE PROCEDURE

**Modified Morrison protocol (Repeated Bead Beating) for DNA extraction from mucosal tissue:**

Reagent and sample prep

1. Thaw sample on ice for 15 minutes (any longer makes it sticky/slimy)
2. Make lysis buffer (can be stored for 1 month)
3. Prefill fresh 2 mL screw cap tubes with 0.4 g sterile zirconia beads (0.2 g of 0.1 mm, 0.2 g of 1 mm)

Cell lysis

1. Add tissue sample to screw-cap tube (use a P20 pipette tip), prefilled with zirconia beads.  
\*Include a "mock" sample
2. Add 300 mL lysis buffer (N.B.: if precipitated, need to heat at 70°C for 15 minutes to fully dissolve).  
Lysis buffer =  
500 mM NaCl,  
50 mM Tris-HCl pH 8.0,  
50 mM EDTA  
4% sodium dodecyl sulphate (SDS).
3. Homogenise in Precelleys – 3 min at 5000 rpm (setting 4)
4. Incubate at 70 °C for 10 minutes
5. Centrifuge for 5 minutes at 16 000 x g (at 4 °C). Transfer supernatant to fresh Eppendorf tube.
6. Add 200 mL lysis buffer then repeat step 3–5. Pool the supernatant

Automated purification

1. Follow instructions on Maxwell/Promega Tissue DNA Purification kit (takes 45 minutes)
2. End product is about 150 – 200 mL of DNA in elution buffer. Transfer to 1.5 mL Eppendorf tube.
3. Centrifuge (1 min @ 5000 x g) DNA in elution buffer to pellet magnetic beads; or place tubes on 1.5 mL magnetic rack and allow beads to collect at side of tube.
4. Aliquot 50 mL of the supernatant into 3 Eppendorfs, label and store in freezer as "gDNA"

Quantification of gDNA<sup>^</sup>:

1. Quantify by nanodrop and record the amount for each sample
2. Visualise on agarose gel (0.8%, 30 minutes, hyperladder 1). Load 2–3 µL.

\*reagents but no tissue, to see if any background bacteria in the DNA extraction/Maxwell procedure (we will sequence and identify this)

^gDNA will be a mixture of 99% human and <1% bacteria. Size is just over 10,000 bp on agarose gel.

### Use of Maxwell for gDNA Extraction from Tissue


- Using samples that have been lysed through bead beating (as per Morrison RBB gDNA extraction protocol above), in approximately 400 µL volume
- Requires Promega Maxwell 16 Tissue DNA Purification kit
- Maxwell SEV mode

Set up of Maxwell machine:

- Must be in SEV mode (see appendix for how to change to SEV).
- Turn on Maxwell. Ensure SEV is shown on home screen.
- Select DNA, then scroll through and select the Tissue DNA protocol.
- Verify that the correct protocol has been chosen, then open door when prompted on the screen.

Preparing Samples:

- Take out required number of cartridges, plungers, elution tubes and elution buffer.
- Label cartridges.
- Place cartridges in preparation rack and peel off seal.
- Pipette your sample (up to 400 µL) into well 1 of the cartridge.
- Place the plunger in well 7 (ridge side) of the cartridge.



- Place the prepared cartridges into the holder in the Maxwell machine (ensure ridged side is clicked in at appropriate end – will only go in one way).
- If you have less than 16 samples, ensure the cartridges are evenly spread across holder starting from the centre.
- Label blue elution tubes corresponding to each of your samples. Place in front of the appropriate cartridge in the elution tube holder.
- Add 300 µL elution buffer to each elution tube.

Running samples:

- Once samples are prepared, press the run button to retract holder. Close the door.
- The run will start automatically. If it is interrupted at any point – your samples will be lost.
- The procedure should take approximately 45 minutes. The time remaining will be displayed on screen.
- Once the run has finished, follow instructions on screen to open door.
- Transfer the contents of the elution tubes into labelled 1.5 mL Eppendorfs.
- Discard the used cartridges and elution tubes.
- Close up machine and set UV to run.

Sample Storage:

- If there is any bead carry over in samples, place tubes on magnetic Eppendorf rack and transfer out supernatant to fresh tubes.
- RNase treat samples if desired.
- Aliquot samples and store at -30 °C.

## **EVALUATION METHOD**

This procedure will be kept updated by the Head of Function Testing, Department of Gastroenterology and Hepatology

## **SUPPORTING DOCUMENTS**

[Measuring bacterial density on biopsy samples using qPCR](#)

## **REFERENCES**

DRAFT